


New Product Release

Replacement parts for
KaVo 6500B Gentlesilence LUX Handpiece




Part Number

HKV8065-A70B2L814

Description

Dynamic Perfection™
Rotor Ceramic



Component Parts

HKV380	Impeller (1)
HKV200	Loading spring washer (1)
HKV205	Loading spring washer (1)
HKV281	Slinger (1)
TSP590	pindle (1)
HKV110	O-Ring (2)
DA02J2L-814	Back Bearing (1)
DA70B2L-801	Front Bearing (1)

**Handpiece
Headquarters**

A  HENRY SCHEIN® COMPANY

For pricing and information regarding
these products, please contact:

Handpiece Headquarters

620 S. Placentia Avenue

Placentia, CA 92870

sales@handpieceheadquarters.com

800-557-7867

714-579-0175

Fax# 714-579-0186