


New Product Release

Dynamic Perfection™ Rotor for
LARES 757MX-ULTRALITE/EURO Handpiece




Part Number

HLR831-R73B3L

Description

Dynamic Perfection Rotor



Component Parts

HLR190	O-ring (1)
HSM166	O-ring (1)
HKV200	Loading Spring (1)
HLR326	Impeller (1)
TSP548-1	Autochuck (1)
DR73B3L	Front Bearing (1)
DR74A3L	Back Bearing (1)

Handpiece Headquarters

2284 N Glassell Street Unit A

Orange CA 92865

Sales@HandpieceHeadquarters.com

Tel: +1-714-685-3070

Fax: +1-714-921-1852

**Handpiece
Headquarters**

A  HENRY SCHEIN® COMPANY

For pricing and information regarding
these products, please contact:

Handpiece Headquarters

620 S. Placentia Avenue

Placentia, CA 92870

sales@handpieceheadquarters.com

800-557-7867

714-579-0175

Fax# 714-579-0186